

AED/Bleed Kit Inspection Form

Inspected By:	Location of Inspection (Building):
Date of Inspection	AED Serial Number:

Use the following maintenance checklist when you check your unit.

Check the following Pass or Fail

Check the following	PASS	FAIL	Comments
Is the green check showing that the unit is ready to use?	<input type="checkbox"/>	<input type="checkbox"/>	
Is the unit clean, undamaged, and free of excessive wear?	<input type="checkbox"/>	<input type="checkbox"/>	
Are there any cracks or loose parts in the housing?	<input type="checkbox"/>	<input type="checkbox"/>	
Verify that electrodes are within their expiration date.	<input type="checkbox"/>	<input type="checkbox"/>	
Verify electrodes are connected to the unit and sealed in their package. Replace if expired.	<input type="checkbox"/>	<input type="checkbox"/>	Exp. Date: